Drempelhulp - Set 1 Montagehandleiding

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
Threshold Ramp – Set 1 Assembly manual

Required tools:


Possibly supplemented by a (screw) drill.

Mounting materials supplied:




Fastening materials, 1 strike plate and 3 laminate feet.

Contents of the box


	part	number	quantity	article number
	Fitting plate, whole 4x3	no. 13	3	8026.100.09
	Fitting plate, half 2x3	no. 14	1	8026.100.10
	Connector 4x1	no. 10	3	8026.100.11
	Connector 2x1	no. 11	1	8026.100.13
	Connector 2x2	no. 12	3	8026.100.12

Optional accessories, not in the box


	Filler cap, grey		14	8025.010.02
	Filler cap, brown		14	8025.010.03
	Anti-skid caps		14	8025.010.01

Assembly of Threshold Ramp

Carefully read the manual before starting the assembly.

Top view

Bottom view

If required, you can adjust the height of the SecuCare Threshold Ramp; you can lower the height of Set 1 by a maximum of 10 mm (in steps of 5 mm). Cutting lines are present in the plastic parts for this.

1 Determine the correct height.



2 Cut to size.

3 If you have had to adjust the height of the sloping part, you will also have to modify the square connector.


4 Place the connectors, square and rectangular, as depicted in the picture below. You can push them in place with your hand, but you can also place the black plastic strike plate on the connector click connection and tap it in place with a hammer. Make sure the surface underneath is solid and flat. Always use the strike plate!

5 The Threshold Ramp can be attached to a slippery surface (e.g. laminate or slippery tiles) using the supplied flat laminate feet. Attach the supplied double-sided tape stickers to the laminate feet.

6 Thoroughly clean / degrease the surface! Attach the laminate feet (place them evenly over the width of the doorstep, at each end and in the middle) and remove the release liner from the double-sided tape.

7 Now place the SecuCare Threshold Ramp in the right location and push it down firmly. Afterwards you can carefully remove the aid and click it back onto the laminate feet.


Attaching the Threshold Ramp to a raised edge/doorstep.

8 Determine where the 3 fastening screws will have to go and cut out the prepared openings.

9 Mark the location on the raised edge/doorstep, pre-drill and fasten the screw and washer. Do not tighten them completely. Now place the SecuCare Threshold Ramp over the screw/washer.

10 The end result:

Disclaimer and warranty conditions – SecuCare Threshold Ramp

Version 3.0 January 2016

Important

- Carefully read the assembly instructions and disclaimer.
- Make sure the surface underneath is solid and flat.
- Be careful not to slip in certain weather conditions. Grit the area in case of freezing rain or do not use the system under such circumstances.
- Make sure that this system does not cause any danger/hindrance (e.g. in a communal walkway).
- The surrounding area must be sufficiently illuminated to ensure that the system is clearly visible, especially in a public place and for people with a visual disability

General provisions

- The Threshold Ramp is installed and used at your own risk.
- The Threshold Ramp is only designed for and suitable for use by persons of an average weight on foot, with a rollator or in a wheelchair, and therefore not for any other application.
- The surface on which the Threshold Ramp is mounted must be flat, stable and solid.
- Never use the system in the vicinity of or in combination with flammable materials. In that case, the system is not suitable for escape routes either.
- The angle of the Threshold Ramp must be as low as possible. If the local situation requires a steep angle, the user(s) must be constantly assisted when using the Threshold Ramp.
- When cutting the Threshold Ramp, adequate protective equipment (safety glasses, ear protection, etc.) must be worn.
- Small parts of the Threshold Ramp must be kept out of reach of children to prevent swallowing, etc.
- When using high heels, it is highly recommended to fit separately available anti-skid caps to the system (art. no. 8025.010.01, 14 pcs). This creates a flat surface (not applicable to set 1).
- After assembly and from time to time afterwards, the Threshold Ramp must be checked to see if it is sturdy enough for normal use in the local situation.
- The area around the doorstep system must be sufficiently illuminated all the time, making the Threshold Ramp clearly visible to others, for whom the Threshold Ramp can be an obstacle after all; for use of the Threshold Ramp in communal walkways, a doorstep system with a sloping side is available.
- Depending on the height of the system, having guide rails (preventing accidents) on both sides of the system is recommended/mandatory. Check your local legislation for this or contact your supplier or Secu for more information.
- Local rules for identifying the Threshold Ramp must be observed. If the Threshold Ramp is slippery, e.g. due to weather conditions, but not limited to this, suitable measures against slipping due to this slipperiness must be taken, whereby the Threshold Ramp can be temporarily removed, for example.
- When in doubt about any aspect of the assembly or use of the Threshold Ramp, the consumer should immediately contact the supplier of the Threshold Ramp from which the consumer purchased it.

Maintenance

- From time to time, the sturdiness of the Threshold Ramp and the individual connections must be checked.
- The Threshold Ramp and the surface underneath must be kept clean at all times. Any type of cleaning product, even an autoclave, can be used to clean the Threshold Ramp.

- Any coarse dirt can be removed with a brush, after which the Threshold Ramp must be rinsed with water.

Warranty

Any warranty claims on the Threshold Ramp require proof of purchase of the Threshold Ramp to be produced, stating the purchase date and product description or product name. Warranty is provided for a period of five years, starting on the date on which the consumer purchased the Threshold Ramp. The warranty only applies to any production errors of the materials used in the Threshold Ramp. A well-founded warranty claim shall only result in the consumer being eligible for Secu supplying a replacement doorstep system free of charge; Secu is not obliged to provide anything else. Secu does not accept any liability whatsoever for any personal injury or any material damage or any consequential loss that occurs during assembly of the Threshold Ramp and/or during use of the Threshold Ramp; any liability to which Secu may be subject despite this is limited to the amount paid by Secu's liability insurer in the relevant case.

The warranty expires immediately if:

- assembly does not meet the assembly instructions supplied with the Threshold Ramp or is performed with materials other than the supplied mounting materials;
- the Threshold Ramp is used incorrectly, improperly or with excessive loads (max. load 850 kg/m²);
- damage is caused to the Threshold Ramp deliberately or through neglect;
- the Threshold Ramp is used on an uneven and/or unstable and/or insufficiently solid surface;
- any damage is caused by a lack of maintenance, by incorrect maintenance or by using the wrong cleaning products (see cleaning recommendation);
- effects of chemical/electrochemical substances and/or effects of electricity occur.

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SecuCare is a trademark of **SECUPRODUCTS BV**