



Rubber Ring

PR20161, PR20618 / 60616, 60618



**Please inspect all parts carefully
for cracks/breaks prior to use.**



Do not over-inflate

USER INSTRUCTIONS


- Features special valve for easy inflation
- Twist open top of valve
- Inflate by mouth or bicycle-style pump to maximum height of 3 inches (7.6 cm)
- Tighten valve to close

CARE INSTRUCTIONS

- Wipe with damp cloth



Maximum
User Wt:



Manufactured For:

Valentine International Ltd., Taiwan R.O.C.

Authorized European Representative:

Medica Surgical Innovations Limited

BB2 4PB, United Kingdom

Made in China

Distributed By:

Able2 UK Ltd.

Moorgate Street, Blackburn, Lancashire, BB2 4PB

Tel: 01254 619000 Fax: 01254 619001

Email: enquiries@able2.eu

Web: www.able2.eu

Rev. 063011